

# Catálogo de productos




## BALANZA WA-30001Y

Rango 3000g x 0,1g



Administración y Ventas: Azcuénaga N° 1627 - Quilmes [ventas@cymed.com.ar](mailto:ventas@cymed.com.ar)  
TE:+54 11 21468039 / 75095121

# HD LCD Display



## High Quality Load Cell



Administración y Ventas: Azcuénaga N° 1627 - Quilmes [ventas@cymed.com.ar](mailto:ventas@cymed.com.ar)  
TE:+54 11 21468039 / 75095121

### PRODUCT DETAILS

Product name	Precision Balance	Description
Item	WA-N	1. Pan: 165×175 mm, stainless iron 304 2. Display: HD LCD 3. Interface: Rs232C 4. Power: AC/DC (Rechargeable Battery is optional) 5. Units: g, kg, oz, ct, lb, GN... 6. More functions: Counting, Percentage, Full Range Tare, Over-Load Protection, Zero Tracking, Automatic Turn-Off, Under Weighing, Low Power Indicator 7. Attachments: english manual, adapter 8. Speed: 2 s 9. Value stable time: <= 2s
Outline Dimension (L×W×H)	192×270×90 mm	
Material breakdown	ABS, STAINLESS IRON, ALUMINIUM ALLOY	
Certificate	CE, ISO, ROHS	
MOQ	12 CARTON BY SHIP	
Country of origin	CHINA	
Port of Loading	SHANGHAI	
H.S. code	9016009000	

### PARAMETER

Model	WA60001N	WA50001Y	WA30001N
Readability	0.1g	0.1g	0.1g
Capacity	0-6000 g	0-5000 g	0-3000 g
Repeatability	<= 0.1 g	<= 0.1 g	<= 0.1 g
Linearity	<= 0.1 g	<= 0.1 g	<= 0.1 g

**CURSOS &  
SEMINARIOS**



[www.cymed.com.ar](http://www.cymed.com.ar)

Administración y Ventas: Azcuénaga N° 1627 - Quilmes [ventas@cymed.com.ar](mailto:ventas@cymed.com.ar)

TE:+54 11 21468039 / 75095121