



BALANZA WA-Y

Rango 2000g x 0,01g, Modelo WA-20002Y

Platillo de 130 x 130 mm

PRODUCT DETAILS

Product name	Precision Balance	Description
Item	WA-Y	1. Pan: 165×175mm (from 2kg to 5kg) Or 130×130mm (from 100g to 2kg), stainless iron 304 2. Display: Big LCD (white back-light) (L×W) 4.3"×1.5" (11×4 cm) 3. Interface: Rs232C 4. Power: AC/DC (Dry Battery) 5. Units: g, kg, oz, ct, lb,... 6. More functions: Counting, Percentage, Full Range Tare, Over-Load Protection, Zero Tracking, Automatic Turn-Off, Under Weighing, Low Power Indicator 7. Attachments: english manual, adapter, calibration weights, under hook (Optional) 8. Speed: 2 s 9. Value stable time: <= 2s
Outline Dimension (L×W×H)	192×270×90 mm	
Material breakdown	ABS, STAINLESS IRON, ALUMINIUM ALLOY	
Certificate	CE, ISO, ROHS	
MOQ	1 CARTON	
Country of origin	CHINA	
Port of Loading	SHANGHAI	
H.S. code	9016009000	

PARAMETER

Model	WA60002Y / WA50002Y	WA40002Y / WA30002Y	WA20002Y / WA10002Y
Readability	0.01g	0.01g	0.01g
Capacity	0-6000 g / 0-5000 g	0-4000 g / 3000 g	0-2000 g / 1000 g
Repeatability	<= 0.02 g	<= 0.02 g	<= 0.01 g
Linearity	<= 0.02 g	<= 0.02 g	<= 0.02 g

Detalles


Dispositivo inteligente	YES	Capacidad	0-5000g
Tamaño	270*192*90mm	Precisión	0.01g
Tipo de pantalla	LCD	Fuente de alimentación	100-240V/50-60HZ
Número de modelo	WA50002Y	Carga nominal	5 kg
Tipo	Báscula de laboratorio	Nombre de la marca	WJEUIP
Origen	CN(Origen)	Certificación	CE,RoHS
Model	WA-Y	Accuracy	0.01g
Capacity	600g-6000g	RS232	Yes
Display	LCD	Calibration Weights	Yes (Capacity 600g to 2000g)



Administración y Ventas: Azcuénaga N° 1627 - Quilmes ventas@cymed.com.ar

TE:+54 11 21468039 / 75095121

Unit transition



Administración y Ventas: Azcuénaga N° 1627 - Quilmes ventas@cymed.com.ar


TE:+54 11 21468039 / 75095121

MORE DETAILS

DETAILS DECIDE
SUCCESS OR FAILURE



Bubble level




Anti-slip feet




RS 232



Load cell protection



INSIDE



Double Pan With Block

Administración y Ventas: Azcuénaga N° 1627 - Quilmes ventas@cymed.com.ar

TE:+54 11 21468039 / 75095121

PRODUCT PARAMETERS



DIMENSION: 192MM×270MM×90MM



● 165×175mm



● 130×130mm



● ϕ 130mm

**CURSOS &
SEMINARIOS**



www.cymed.com.ar

Administración y Ventas: Azcuénaga N° 1627 - Quilmes ventas@cymed.com.ar

TE:+54 11 21468039 / 75095121